

HOLTER SYSTEM

EC-3H/ABP COMBINED HOLTER SYSTEM (BLUETOOTH)


- **FULL-DISCLOSURE AMBULATORY ECG AND BP RECORDS**
- 3 INDEPENDENT ECG CHANNELS
- OSCILLOMETRIC BP MEASUREMENT
- **PACEMAKER** CHANNEL AND ANALYSIS
- PATIENT EVENT BUTTON
- RECHARGEABLE BATTERIES FOR ECONOMIC OPERATION
- USER-FRIENDLY SOFTWARE
- FAST QRS CLASSIFICATION AND RHYTHM ANALYSIS
- **HRV, ST, PQ, QT, PM, ABP AND ATRIAL FIBRILLATION, FLUTTER ANALYSIS, BBB AND AV BLOCK DETECTION**
- COLOR-CODED GRAPHS AND TABLES
- VARIOUS EDITING AND PRINTING FUNCTIONS




EC-3H/ABP



The EC-3H/ABP combines a 3-channel ECG Holter and an ambulatory blood pressure monitor. The QRS and PM can be detected with high precision, while template and rhythm analysis can be interactively modified by the user. This user-friendly software provides high-quality time- and frequency-based HRV, PQ, QT, ST, PM, Atrial Fibrillation, Flutter analysis, BBB and AV block detection and ABP results. The transparency and the interpretation of the results are aided by color-coded graphs and tables. The software offers a wide range of languages to choose from.



FULL VIEW + BP



ST 3D